



Biosolve specialist of high purity solvents for LC-MS, recognized and recommended by the manufacturers of chromatograph, offer you:

ACETONITRILE  
**ULC/MS**  
**CC/SFC**  
 CODE #: 012041

New assay:  
**99.99%**

- Carton box 4 x 2.5 or 6 x 1 litres
- Glass bottle
- Packed under nitrogen with new filtration through 0.03 µm



To always raise the bar

**higher!**

# ACETONITRILE

## ULC/MS – CC/SFC

CODE #: 012041

Test Description	Specifications
Appearance	Clear colorless liquid
Color (APHA, Pt-Co)	0-5
Assay (GC, on anhydrous basis)	99.99-100 %
Residue after evaporation	0-0.0001 %w/w
Water (KF)	0-0.01 %w/w
Acidity (as Acetic acid)	0-0.001 %
Alkalinity (as Ammonia)	0-0.0001 %
MS-ESI+ (as Reserpine)	0-2 ppb
H.PEAK by PDAD 210-400nm	0-0.001 AU
Grad. elution H.PEAK at 210nm	0-0.001 AU
Grad. elution drift at 210nm	0-0.006 AU
Grad. elution H.PEAK at 254nm	0-0.0003 AU
Grad. elution drift at 254nm	0-0.002 AU
F254nm (as Quinine)	0-0.3 ppb
F365nm (as Quinine)	0-0.3 ppb
T191nm	30-100 %
T195nm	85-100 %
T200nm	97-100 %
T215nm	98-100 %
T>230nm	99-100 %

- Carton box 4 x 2.5 or 6 x 1 litres
- Glass bottle
- Packed under nitrogen with new filtration through 0.03 µm

Contact us  
with your annual  
needs **to have**  
**the best price.**

CIENTISOL, SLU  
T. 981 936 338  
info@cientisol.com · www.cientisol.com



Because you  
need **Quality**

